

Accessories for Electrophoresis

Gradient Mixers, Gel Dryer



Gradient Mixer

- for preparation of polyacrylamide gel gradients or centrifugation gradients
- made of robust acrylic block
- in 5 different sizes available

Description	Order No.
Gradient mixer, 15 ml	GM15
Gradient mixer, 25 ml	GM25
Gradient mixer, 50 ml	GM50
Gradient mixer, 100 ml	GM100
Gradient mixer, 500 ml	GM500



Gel Dryer

- even heat distribution via an 800 W heating membrane
- dries sequencing gels up to 40 x 50 cm in less than 30 minutes
- accurate temperature control up to 90 °C
- clear silicone rubber overlay sheet
- 5 hour timer for heating
- connection to vacuum pump possible
- 5 hour timer for vacuum pump

Description	Order No.
45 x 34 cm Gel Dryer including stainless steel screen, mylar sheet, porous polyethylene sheet, clear rubber overlay sheet	GD-4534
50 x 40 cm Gel Dryer including stainless steel screen, mylar sheet, porous polyethylene sheet, clear rubber overlay sheet	GD-5040