

Accessories for Automated Sequencing Units

Suitable for LI-COR® sequencing systems

Gelco® Rectangle combs

Order No.	Description	Quantity
60-11-107	34well rectangle comb, mylar, 0.20 mm thick	1 piece
60-11-147	34well rectangle comb, mylar, 0.25 mm thick	1 piece
60-11-148	34well rectangle comb, mylar, 0.40 mm thick	1 piece
60-11-108	48well rectangle comb, mylar, 0.20 mm thick	1 piece
60-11-149	48well rectangle comb, mylar, 0.25 mm thick	1 piece
60-11-150	48well rectangle comb, mylar, 0.40 mm thick	1 piece
60-11-156	50well rectangle comb, mylar, 0.20 mm thick	1 piece
60-11-157	50well rectangle comb, mylar, 0.25 mm thick	1 piece
60-11-158	50well rectangle comb, mylar, 0.40 mm thick	1 piece
60-11-109	64well rectangle comb, mylar, 0.20 mm thick	1 piece
60-11-151	64well rectangle comb, mylar, 0.25 mm thick	1 piece
60-11-152	64well rectangle comb, mylar, 0.40 mm thick	1 piece
60-11-159	66well rectangle comb, mylar, 0.20 mm thick	1 piece
60-11-160	66well rectangle comb, mylar, 0.25 mm thick	1 piece
60-11-161	66well rectangle comb, mylar, 0.40 mm thick	1 piece
60-11-153	72well rectangle comb, mylar, 0.20 mm thick	1 piece
60-11-154	72well rectangle comb, mylar, 0.25 mm thick	1 piece
60-11-155	72well rectangle comb, mylar, 0.40 mm thick	1 piece

Gelco® Sharkgum combs

Order No.	Description	Quantity
60-11-129	32well sharkgum, 0.20 mm thick	1 piece
60-11-130	36well sharkgum, 0.20 mm thick	1 piece
60-11-131	48well sharkgum, 0.20 mm thick	1 piece
60-11-132	64well sharkgum, 0.20 mm thick	1 piece
60-11-133	72well sharkgum, 0.20 mm thick	1 piece
60-11-126	96well sharkgum, 0.20 mm thick	1 piece

Gelco® Casting combs

Order No.	Description	Quantity
60-11-128	Casting comb, 0.25 mm thick	1 piece

LI-COR® is a registered trademark of LI-COR Inc. USA.

Note: About special combs

Sharkgum Combs

- have square bottom teeth but are actually used as a typical sharkstooth comb