

Accessories for Protein Electrophoresis

Suitable for BIOMETRA™ Minigel and Minigel Twin

Gelco® Glass plates

Order No.	Description	Quantity
50-30-102	Glass plate, 10.5 x 9.8 cm with spacer 1.0 mm thick	1 piece
50-30-103	Glass plate, 10.5 x 9.8 cm, notched	1 piece

Gelco® Combs

Order No.	Description	Quantity
50-31-103	5well teflon comb, 1.00 mm thick	1 piece
50-31-106	10well teflon comb, 1.00 mm thick	1 piece
50-31-109	15well teflon comb, 1.00 mm thick	1 piece
50-31-112	20well teflon comb, 1.00 mm thick	1 piece

Suitable for BIOMETRA™ Multigel

Gelco® Glass plates

Order No.	Description	Quantity
50-30-112	Glass plate, 14.8 x 9.3 cm with spacer 1.0 mm thick	1 piece
50-30-113	Glass plate, 14.8 x 9.3 cm, notched	1 piece

Suitable for BIOMETRA™ Multigel-Long

Gelco® Glass plates

Order No.	Description	Quantity
50-30-122	Glass plate, 14.8 x 14.3 cm with spacer 1.0 mm thick	1 piece
50-30-123	Glass plate, 14.8 x 14.3 cm, notched	1 piece

Suitable for BIOMETRA™ Maxigel

Gelco® Glass plates

Order No.	Description	Quantity
50-30-131	Glass plate, 20.8 x 20.2 cm with spacer 1 mm thick	1 piece
50-30-132	Glass plate, 20.8 x 20.2 cm with spacer 2 mm thick	1 piece
50-30-133	Glass plate, 20.8 x 20.2 cm, notched	1 piece

Biometra® is a registered trademark of Biometra biomedizinische Analytik GmbH.