

# Pipette Tips

## Standard tips

biostep universal tips are engineered for a wide variety of pipettes and are ideal for laboratories with multiple brands of pipettes – even multi channel pipettes. The universal tips feature the universal grip design, using flexible materials with a series of bands and supports producing a positive seal. biostep tips have thin walls and a small surface area at the tip, greatly minimizing the possibility of the tip retaining template droplets.



In addition, the small surface area reduces the touch-off angle between the outer tip wall and the receiving sample tubes. This reduced touch-off angle virtually eliminates the potential of capillary retention between the tip orifice and tubes wall. All tips are provided with calibration lines at specific volumes, acting as visual checks to minimize any risk of errors during sample aspiration.

### 0.1 to 10 µl, natural

---

pipet4u digital micro pipettes  
Gilson Pipetman P-2/P-10  
Eppendorf Reference/Research  
Labsystems Finn pipettes  
Biohit Proline pipettes  
Nichiryo pipettes  
Oxford Benchmate pipettes  
Rainin EDP-2 pipettes  
HTL pipettes



### 0.5 to 10 µl, natural

---

pipet4u digital micro pipettes  
Gilson pipettes  
Eppendorf Reference/Research  
Nichiryo pipettes  
Oxford Benchmate pipettes  
Socorex pipettes  
Labsystems Finn pipettes  
HTL pipettes



### 200 µl, natural and yellow

---

pipet4u digital micro pipettes  
Eppendorf Reference/Research  
Gilson pipettes  
Brand pipettes  
Labsystems Finn pipettes  
Biohit Proline pipettes  
Socorex pipettes  
Nichiryo pipettes  
MLA pipettes  
Costar pipettes  
Oxford Benchmate pipettes  
Rainin EDP-2 pipettes  
HTL pipettes



### 300 µl, natural

---

pipet4u digital micro pipettes  
Eppendorf Reference/Research  
Gilson pipettes  
Brand pipettes  
Labsystems Finn pipettes  
Biohit Proline pipettes  
Oxford Benchmate pipettes

# Pipette Tips

## Standard tips

### 100 to 1000 $\mu$ l, natural and blue

pipet4u digitale micro pipettes  
Eppendorf Reference/Research  
Gilson pipettes  
Brand pipettes  
Labsystems Finn pipettes  
Biohit Proline pipettes  
Socorex pipettes  
Nichiryo pipettes  
MLA pipettes  
Costar pipettes  
Oxford Benchmate pipettes  
HTL-pipettes  
Rainin EDP-2 pipettes



### Tips in bags, non-sterile

Description	Colour	Quantity	Order No.
0.1 - 10 $\mu$ l	natural	1 bag with 1000 pieces	01-10-101
0.5 - 10 $\mu$ l	natural	1 bag with 1000 pieces	01-10-111
200 $\mu$ l	natural	1 bag with 1000 pieces	01-10-121
200 $\mu$ l	yellow	1 bag with 1000 pieces	01-10-131
300 $\mu$ l	natural	1 bag with 1000 pieces	01-10-151
1000 $\mu$ l, long shank	natural	2 bags with 500 pieces	01-10-161
1000 $\mu$ l, long shank	blue	2 bags with 500 pieces	01-10-171
1000 $\mu$ l, long shank	blue	1 bag with 1000 pieces	01-10-191

### Tips in racks, non-sterile

Description	Colour	Quantity	Order No.
0.1 - 10 $\mu$ l	natural	10 racks with 96 pieces	01-10-102
0.5 - 10 $\mu$ l	natural	10 racks with 96 pieces	01-10-112
200 $\mu$ l	natural	10 racks with 96 pieces	01-10-122
200 $\mu$ l	yellow	10 racks with 96 pieces	01-10-132
300 $\mu$ l	natural	10 racks with 96 pieces	01-10-152
1000 $\mu$ l, long shank	natural	10 racks with 96 pieces	01-10-162
1000 $\mu$ l, long shank	blue	10 racks with 96 pieces	01-10-172
1000 $\mu$ l, long shank	blue	10 racks with 60 pieces	01-10-192

### Tips in racks, sterile

Description	Colour	Quantity	Order No.
0.1 - 10 $\mu$ l	natural	10 racks with 96 pieces	01-10-202
0.5 - 10 $\mu$ l	natural	10 racks with 96 pieces	01-10-212
200 $\mu$ l	natural	10 racks with 96 pieces	01-10-222
300 $\mu$ l	natural	10 racks with 96 pieces	01-10-228
1000 $\mu$ l, long shank	natural	10 racks with 96 pieces	01-10-232
1000 $\mu$ l, short shank	blue	10 racks with 60 pieces	01-10-242