





www.biostep.de
www.techne.de

Consumables for every lab ... an all-round solution



- good quality for low budget
- short delivery times
- products made of PP are guaranteed free from DNase, RNase, pyrogen and metal
- within Germany: no freight charges when ordering for more than €100.00

- PCR tubes and PCR plates
- reaction tubes
- tips and filter tips
- micro pipettes
- adhesive foils
- strips
- labels
- and much more ...



Your Specialist for PCR, Electrophoresis, Bio-Imaging & Radioanalytics

PIPETTING

- tips and filter tips for each laboratory
- universal fit
- for most commonly used pipettes on the market
- made of polypropylene (PP)
- produced under clean room environments
- free from RNase, DNase and pyrogen
- autoclavable at 121 °C, 0.1 Mpa, 20 min



Tips in bags, non-sterile

Description	Quantity	Order No.	Price in Euro
0.1 - 10 µl natural	1 bag of 1000 pieces	01-10-101	10.00
0.5 - 10 (20) µl natural	1 bag of 1000 pieces	01-10-111	10.00
200 µl natural	1 bag of 1000 pieces	01-10-121	7.00
200 µl yellow	1 bag of 1000 pieces	01-10-131	7.00
300 µl natural	1 bag of 1000 pieces	01-10-151	12.00
1000 µl natural, long	1 bag of 1000 pieces	01-10-161	10.00
1000 µl blue, long	1 bag of 1000 pieces	01-10-171	10.00
1000 µl blue, short	1 bag of 1000 pieces	01-10-191	10.00

Filter tips in racks, sterile

Description	Quantity	Order No.	Price in Euro
0.1 - 10 µl natural	10 racks of 96 pieces	01-11-202	50.00
0.5 - 10 µl natural	10 racks of 96 pieces	01-11-212	50.00
200 µl natural	10 racks of 96 pieces	01-11-262	50.00
230 µl natural	10 racks of 96 pieces	01-11-272	55.00
1000 µl natural, long	10 racks of 96 pieces	01-11-282	55.00
1000 µl blue, short	10 racks of 60 pieces	01-11-292	50.00



PIPETTES

- high accuracy and precision
- approved push-button mechanism
- easy-to-read display with large numbers
- convenient one-hand operation with ergonomic soft grip
- 2-step plunger operation, blow-out system
- simple maintenance and calibration




Micro pipettes, digital volume display, fixed and variable volumes

Description	Volume	Quantity	Order No.	Price in Euro
Micro pipette, fixed volume	10 µl	1 piece	PIPET-F31	49.00
Micro pipette, fixed volume	100 µl	1 piece	PIPET-F35	49.00
Micro pipette, fixed volume	1000 µl	1 piece	PIPET-F39	49.00
Micro pipette, variable volume	0.5 to 10 µl	1 piece	PIPET-D10	89.00
Micro pipette, variable volume	10 to 100 µl	1 piece	PIPET-D13	89.00
Micro pipette, variable volume	100 to 1000 µl	1 piece	PIPET-D16	89.00
8 channel pipette, variable volume	0.5 to 10 µl	1 piece	PIPET-M60	349.00
8 channel pipette, variable volume	5 to 50 µl	1 piece	PIPET-M61	349.00
8 channel pipette, variable volume	30 to 300 µl	1 piece	PIPET-M62	349.00



PCR

- **PCR tubes and PCR plates** for the use in thermal cyclers of the most well-known manufactures on the market
- colour: natural, on request also in other colours available e.g. blue, green, red, yellow...
- free from RNase, DNase, pyrogen and metal
- made of polypropylene (PP)
- perfectly suitable for all commonly used thermal cyclers
- extremely thin walled for a uniform and fast temperature transfer



PCR tubes 0.2 and 0.5 ml with caps

Description	Quantity	Order No.	Price in Euro
0.2 ml tubes flat cap	2 x 500 pieces	01-01-101	37.00
0.2 ml tubes domed cap	2 x 500 pieces	01-01-104	37.00
Rack for 96 x 0.2 ml tubes, yellow	1 piece	01-01-510	3.95
0.5 ml tubes flat cap	2 x 500 pieces	01-01-201	32.00
0.5 ml tubes domed cap	2 x 500 pieces	01-01-204	32.00
Rack for 60 x 0.5 ml PCR tubes, orange	1 piece	01-01-520	3.95

PCR tube strips, 8 connected 0.2 m tubes with separate cap strips

Description	Quantity	Order No.	Price in Euro
0.2 ml strips of 8 tubes with cap strips, flat	250 pieces ea.	01-01-304	169.00
0.2 ml strips of 8 tubes with cap strips, flat, LP**	250 pieces ea.	01-01-403	165.00
0.2 ml strips of 8 tubes with cap strips, domed	250 pieces ea.	01-01-302	172.00
0.2 ml strips of 8 tubes with cap strips, domed, LP**	250 pieces ea.	01-01-402	168.00

**LP = Low Profile: reduced height of the tubes for reduction of the dead volume

PCR plates, 48/96well, alphanumeric label

Description	Quantity	Order No.	Price in Euro
Thermo-Fast® 48well plate, semi-skirted	50 pieces	01-03-050	99.00
Thermo-Fast® 96well plate, non-skirted	25 pieces	01-03-101	76.00
Thermo-Fast® 96well plate, semi-skirted	25 pieces	01-03-111	76.00
Thermo-Fast® 96well plate, semi-skirted with cross web	25 pieces	01-03-121	76.00
Thermo-Fast® 96well plate, semi-skirted, detection	25 pieces	01-03-131	86.00
Thermo-Fast® 96well plate, non-skirted, LP**	25 pieces	01-03-201	76.00
Thermo-Fast® 96well plate, skirted, LP**	25 pieces	01-03-211	76.00

**LP = Low Profile: reduced height of the tubes for reduction of the dead volume


Cap strips and self-adhesive foils for sealing PCR plates

Description	Quantity	Order No.	Price in Euro
Cap strips of 8, domed	250 pieces	01-01-501	42.00
Cap strips of 8, flat	250 pieces	01-01-502	42.00
Adhesive aluminium foil, from -70 °C to 150 °C	100 sheets	01-05-101	70.00
Adhesive transparent foil, from -20 °C to 120 °C	100 sheets	01-05-102	86.00

suitable for QPCR

Cap strips of 8, ultra clear, flat	120 pieces	01-01-602	46.00
Adhesive transparent foil, from -40 °C to 120 °C	50 sheets	01-05-201	82.00

The listed trademarks are property of the respective companies.





LABELING

- pre-cut, peel-off labels for 0.5 to 2.0 ml tubes
- temperature resistant from -40 °C to +120 °C
- suitable et al. for the use in boiling water baths, for centrifugation and steam-sterilisation
- 1000 pieces per roll
- 23 x 13 mm labels are suitable for 0.5 ml tubes
- 32 x 13 mm labels are suitable for 1.5 ml and 2.0 ml tubes

Labels, tearproof

Description	Quantity	Order No.	Price in Euro
Tough-Tag™ label 23 x 13 mm, red	1 roll of 1000 pieces	TTSR1000	45.00
Tough-Tag™ label 23 x 13 mm, yellow	1 roll of 1000 pieces	TTSY1000	45.00
Tough-Tag™ label 23 x 13 mm, blue	1 roll of 1000 pieces	TTSB1000	45.00
Tough-Tag™ label 23 x 13 mm, green	1 roll of 1000 pieces	TTSG1000	45.00
Tough-Tag™ label 32 x 13 mm, red	1 roll of 1000 pieces	TTLR1000	45,00
Tough-Tag™ label 32 x 13 mm, yellow	1 roll of 1000 pieces	TTLY1000	45.00
Tough-Tag™ label 32 x 13 mm, blue	1 roll of 1000 pieces	TTLB1000	45.00
Tough-Tag™ label 32 x 13 mm, green	1 roll of 1000 pieces	TTLG1000	45.00

The listed trademarks are property of the respective companies.



CENTRIFUGATION

- tubes made of polypropylen (PP) with graduation
- free from RNase, DNase, pyrogen and heavy metal
- with marking field on the side and writeable cap
- max. RCF for
micro centrifuge tubes: 18.000 x g (fixed-angle rotor)
centrifuge tubes: 9.500 x g (fixed-angle rotor)
- autoclavable without cap

Micro centrifuge tubes, 1.5 and 2.0 ml with cap, non-sterile

Description	Quantity	Order No.	Price in Euro
1.5 ml tubes, with cap, clear	1 pack of 1000 pieces	01-08-103	15.00
2.0 ml tubes, with cap, clear	1 pack of 1000 pieces	01-08-201	19.00

Centrifuge tubes, 15 and 50 ml with blue screw cap, sterile

Description	Quantity	Order No.	Price in Euro
15 ml centrifuge tubes	1 pack of 500 pieces	01-08-302	71.00
50 ml centrifuge tubes	1 pack of 500 pieces	01-08-402	79.00
50 ml centrifuge tubes, with skirt	1 pack of 450 pieces	01-08-403	91.00



FREEZING

- free from RNase, DNase, DNA and pyrogen
- temperature resistant from -90 °C to +130 °C
- with white marking field

Cryo tubes, 1.8 ml with cap and silicon „o-ring mounted“

Description	Quantity	Order No.	Price in Euro
1.8 ml cryo tubes	1 pack of 500 pieces	01-08-701	45.00
1.8 ml cryo tubes, sterile	1 pack of 500 pieces	01-08-702	62.00
Rack, clear/magenta	1 pack of 5 pieces	01-08-750	15.00

Please request your sample - free of charge:

- Tips and filter tips
- PCR tubes
- Cryo tubes



www.biostep.de • www.techne.de

biostep® GmbH
 Meinersdorfer Straße 47A
 09387 Jahnsdorf, Germany
 Tel. +49(0) 3721-39 05-0
 Fax +49(0) 3721-39 05-28
 E-Mail: vbm@biostep.de